Exhibit A



Hamilton Medical, Inc. | 201 Edison Way, Unit A | Reno | NV 89502-2305 | USA phone +1 800 426 6331 | email info@hamilton-medical.com | website hamilton.medical.com

September 17th, 2024

Subject: HAMILTON-T1 transport ventilator

To Whom It May Concern:

Hamilton Medical is the sole source manufacturer and distributor of the HAMILTON-T1 transport ventilator. Hamilton Medical sells directly to the end-user facility and as such, there are no agents or dealers authorized to represent this product. No other division of Hamilton Medical, nor any other company, makes similar or competing products.

Key characteristics:

- Dynamic Lung
 - o Provides direct visualization of real-time changes in pulmonary mechanics.
- Closed Loop Ventilation with Adaptive Support Ventilation[®] (ASV)
 - Provides automated selection and modification of respiratory rate (RR), tidal volume (Vt), and spontaneous breathing resulting in:
 - Lower respiratory system driving pressure
 - Lower inspiratory transpulmonary pressure delivery
 - Lung protective ventilation
- Airway Pressure Release Ventilation (APRV)/DuoPap modes
 - Improved oxygenation
 - Preservation of spontaneous breathing
- Volumetric Capnography (VCO2)
 - Hamilton Medical Volumetric CO2 technology provides direct visualization of the heart/lung interaction to guide optimal ventilation strategy based on patient physiology (other companies only measure and display End Tidal CO2)
- IntelliTrig®
 - Standard feature associated with reducing Patient/Ventilator trigger asynchrony, ultimately resulting in reduction in Asynchrony Index (AI)

